

Go Beyond the Scale

Take an InBody Test Today!

In just 30 seconds, the InBody Result Sheet will show your body fat, muscle, and body water, giving you a better understanding of your weight and health.

The InBody 300 is a sleek, white and silver body composition analyzer. It features a central vertical panel with the 'InBody' logo and two black electrode pads at the bottom. A red arrow points from the left side of the device to a list of six features, each with a red circular icon. Another red arrow points from the bottom right of the device to a sample 'InBody Score' result sheet.

- Muscle Mass
- Body Water
- Segmental Fat
- Percent Body Fat
- Visceral Fat Level
- Basal Metabolic Rate
- History
- Phase Angle

InBody Score 67

Body Composition Analysis

Parameter	Value	Unit
Weight	64.6	kg
Extracellular Water	24.3	(kg)
Dry Lean Mass	21.6	(kg)
Body Fat Mass	21.6	(kg)
Muscle Fat Analysis	47.8	(kg)

Obesity Analysis

Parameter	Value	Unit
BMI	22.0	kg/m ²
Waist Circumference	80.0	cm
PPG	100	ml/min

Segmental Lean Analysis

Region	Right Arm	Left Arm	Trunk	Right Leg	Left Leg
Weight	1.9	1.9	1.9	1.9	1.9
Extracellular Water	0.5	0.5	0.5	0.5	0.5
Dry Lean Mass	1.4	1.4	1.4	1.4	1.4
Body Fat Mass	0.5	0.5	0.5	0.5	0.5

ECW/TBW-Phase Angle

Region	Right Arm	Left Arm	Trunk	Right Leg	Left Leg
ECW/TBW	0.290	0.290	0.290	0.290	0.290
Phase Angle	40.0	40.0	40.0	40.0	40.0

Body Composition History

Parameter	Current	1 Month Ago	3 Months Ago	6 Months Ago	1 Year Ago
Weight	64.6	64.6	64.6	64.6	64.6
BMI	22.0	22.0	22.0	22.0	22.0
PPG	100	100	100	100	100
ECW/TBW	0.290	0.290	0.290	0.290	0.290

Results Interpretation QR Code

Scan the QR Code to see your results interpreted by a healthcare professional.

Impedance

Impedance (Ω) vs. Region (RA, LA, TR, RL, LL, LL)

RA: 12.0, LA: 12.0, TR: 12.0, RL: 12.0, LL: 12.0, LL: 12.0

INBODY 300. © 2020 INBODY USA. INBODY is a registered trademark of INBODY USA, Inc. All rights reserved.